

ROOTT implant library for Exocad software



We recommend following these rules :

- Do not mix different products between implant categories. The only exception is the ROOTT R library. In which ROOTT S and ROOTT M scan posts with M1/MS1 abutments could be used with ROOTT R library.
- ROOTT R, ROOTT M, ROOTT S, ROOTT P implant systems intended for two-piece implants, ROOTT C, ROOTT CS, ROOTT B, ROOTT BS – for one piece.
- Two-piece implants work in principle following implant + together with an abutment, and one-piece implants are based on the telescopic abutment.
- Library name are created from terms: Manufacturer+ Implant + abutment, material.
- Be sure about correct digital and scan post position—it guarantee workflow accuracy(how to detect correct identification).
- Products can only be used as long as the product is in pristine/precise condition: the surface is not scratched, reshaped, deformed (any change in the shape of the scanning body may adversely affect the accuracy of the scanning process and the subsequent accuracy of the workflow and result).

The details included in Exocad library

- ROOTT B; ROOTT BS; ROOTT C; ROOTT CS (External platform telescopic abutments: PCE0 – PCE3, PCES1–PCES3, PCEXS1–PCEXS2, TCE0 – TCE3, TCES1–TCES3, TCEXS0–TCEXS2)
- ROOTT M (PCOM, digital detail for framework)
- ROOTT P (PCOM, digital detail for framework)
- ROOTT S (PCOMS, digital detail for framework)
- ROOTT R (PCO1 – PCO3, POC1S – PCO3S, PCO, CRE, PMAB)
- Scan post (SPCO, SPCOIO, SPCOM, SPCOMIO, SPCOMS, SPCOMIOS, TRA, HE, TOEA, TOE, External platform)
- Digital analogs (AND, ANMD, ANMSD, ANED).

DISCLAIMER OF WARRANTIES.

ALL SERVICES AND TRATE MATERIALS ARE PROVIDED “AS IS” AND “AS AVAILABLE” WITHOUT WARRANTY OF ANY KIND, AND TRATE HEREBY DISCLAIMS ALL WARRANTIES, WHETHER EXPRESS, IMPLIED, STATUTORY, OR OTHER. TRATE SPECIFICALLY DISCLAIMS ALL IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, TITLE, AND NON-INFRINGEMENT, AND ALL WARRANTIES ARISING FROM THE COURSE OF DEALING, USAGE, OR TRADE PRACTICE.

WITHOUT LIMITING THE FOREGOING, TRATE MAKES NO WARRANTY OF ANY KIND THAT THE SERVICES, TRATE MATERIALS, OR THE INFORMATION CONTAINED IN THE SERVICES WILL MEET THE USER’S OR ANY OTHER PERSON’S REQUIREMENTS, OPERATE WITHOUT INTERRUPTION, ACHIEVE ANY INTENDED RESULT, BE COMPATIBLE OR WORK WITH ANY SOFTWARE, SYSTEM, OR OTHER SERVICES, OR BE SECURE, ACCURATE, TIMELY, RELIABLE, RELEVANT, COMPLETE, FREE OF HARMFUL CODE, OR ERROR-FREE.

ROOTT does not define the maximum number of uses appropriate for reusable san posts. The useful life of these devices depends on a number of factors including the methods and of each uses and the handling between uses. Products should not be used if these defects are visible:

- Corrosion, rusting;
- Pitting, discoloration


For precise and quality scanning, you must take these factors into account:

- Accurate scanning body geometry and position
- Surface and shape of the scanning body
- Number of frames made per second
- The ability of a scanning machine to calculate similar surfaces and join them together



Scan Posts

Scan posts indicate the exact position of the implant in the jaw. During the scanning process, the information about the position is transferred into digital format. ROOTT library has intraoral and extraoral scan posts. The difference between them is size and shape, which provide better performance for a particular workflow. The scheme below shows suitable scan posts for a specific type of implant. Working on with ROOTT R implants, scan posts (SPCOM, SPCOMIO, SPCOMS, SPCOMIOS) can be used by assembling them with abutments M1 or MS1. Another advantage of the digital ROOTT library is that transfer, and telescopic abutments of ROOTT C, CS, B, BS implants also can be used as scan posts. These mentioned possibilities create broader applicability of the products.




Precise detection of scan-bodies

Orientation of scan data shall be chosen correctly, whereas it determines how precisely the program detects the scan body.


Reason for undetected scan body:

- Scan data is not precise;
- The scanner is not calibrated;
- Scan body is damaged or defected;
- Scan body is screwed incorrectly;
- Poor Intraoral scanner's scanning quality


Inaccurate detection results



1.SPCO scan post. Good quality of scan and detection.




2.SPCO scan post. The medium level of scanning quality, perfect quality of detection.




3. Apt SPCOIO scan post.

Presices detection results



1.SPCOIO scan post. Damaged surfaces



2. CRE.06 abutment and scan post. Damaged surfaces



3. TRA scan post, Inaccurate scan.

Cement gap

An empty gap between a titanium base and an inner crown surface is called a cement gap. A dimension in microns (μm) indicates an offset from the titanium base.

Cement gap size results depend on the equipment. To get the best results, it is recommended:

- to find the most suitable cement gap option for your equipment;
- to check regularly the mill in the machine;
- to comply with the requirements of the equipment manufacturers;
- to check if the setting of the prosthesis and manufacturing equipment match (see picture below).




ROOTT digital library is created for 3 cement gap categories:

- Zirconia - 0.025, 0.035, 0.050 μm ;
- Plastic PEEK or PMMA - 0.075, 0.090, 0.110 μm .
- Metal - 0.015, 0.025 μm .

If there is a need to create from PEEK/PMMA with a narrow cement gap, choose Zro section.

Metal cement gaps is adapted only for ROOTT M and S Multi Unit and External platform - Rot peek pmma/zro.




Screw channel


There are two options to create an aesthetical prosthesis:

- Standard screw channel that matches the diameters of the screw and the screwdriver.
- Narrow screw channel that matches only the screwdriver. Choose this option if there is less occlusal surface, or a more aesthetical result is desired.


Screw channel sizes are available for ROOTT R, ROOTT M, ROOTT P, and ROOTT S.



Angled screw channel position can be determined by scan post. The front plane of the scan post corresponds to the angle direction of the screw channel. Therefore, we recommend turning the scan post to lingual surface or palatal surface direction before scanning in order to achieve the best results.



The front plane of scan posts



Correct scan post ROOTT R direction highlighted green, incorrect - red



Abutments

ROOTT R digital abutments are suitable for ordinary and complicated clinical situations. There is a wide range of options for multi-unit and single crown cases. In this instruction You could see all abutment which are included in to digital library.

CRE is a multi-functional part that is made of the same material as an implant and abutment. It is applicable as an abutment for immediate loading, open/close tray transfer, carrier for implant insertion, or healing abutment.

For multi-unit cases, use M1+PCOM or M1 Multi-Unit and for a single crown, choose from a PCO titanium base.

Pre-milled abutment PMAB is a customizable, one-piece abutment for a single crown metal or plastic PEEK framework.

ROOTT R abutment characteristics

- Angled access for tunnel from 0° to 20°;
- For bridges (PCOR, MS1, M1, CRE,) and single crowns (PCO1-PCO3, PMAB);
- Variety of gingiva part height PCO1-PCO3S;
- Variety of titanium bases height PCO and PCOS.

PMAB abutment characteristics

- Approved for use with a MEDENTIKA PreFace® Abutment Blank Holders;
- Provides unlimited possibilities to create high precision one-piece customized titanium abutment with an in house milling machine;
- Ideal adjustment for shape, emergence profile, esthetic properties - are available for frequent situations.



ROOTT R

Scan posts



Scan posts indicate the exact position of the implant in the jaw. During the scanning process, the information about the position is transferred into digital format. According to scan post position, height and direction, scan post is converted to abutment and sets the precise location of analogue. Get yourself acquainted with all types of ROOTT R scan posts, developed for the effective and precise workflow of the Odontology specialist (intraoral) and dental technicians (extraoral).

	 SPCO	 SPCOIO	 M1		 MS1	
			 SPCOM	 SPCOMIO	 SPCOMS	 SPCOMIOS
Long body	✓		✓		✓	
Short body		✓		✓		✓
Easily scannable	✓	✓	✓	✓	✓	✓
Positionating	✓	✓				
Rotational			✓	✓	✓	✓
Using with all ROOTT R libraries	✓	✓				
Using with ROOTT R_M1/MS1 libraries			✓	✓	✓	✓

ROOTT R Two-Piece Implant > Internal cone > Titanium base



Same library PCO1/ PCOR- few differences


- Same gingiva and ti base height;
- PCOR without hexagon has been specifically designed for bridges to solve slight divergence between more than one implant. PCO1 was specifically created for single units to help prevent rotation and create stability for the screw.
- More usage options with combining these parts together.
- PCO1-3 Titanium bases with flat engaging spot on the body of the base near the platform.



Non- Enganging - PCOR



Enganging- PCO1, PCO2, PCO3 and others


































































More usage options with combining these parts together, options such as:

- Use PCO1
- When You making single crown or bridge.
- Use PCOR when You making more than one unit prosthesis.
- Combine PCO1 or any other PCOx detail together with PCOR and achieve burden-free prosthesis consturction fixation. If both interfaces had the hex in some cases it would be impossible to get the ti-base to seat on both implants.
- Another difference is the flat engaging spot on the body of the base near the platform.



ROOTT library file system for Exocad software


1. Library folders start with the Manufacturer's name- Trate.
2. ROOTT System consists of:
 - ROOTT External platform - libraries for ROOTT C/ ROOTTCS/ ROOTT B/ ROOTT BS.
 - ROOTT R- two-piece implant system library
 - ROOTT MP/S- Multi-unit implant system library
 - ROOTT K - Coming soon*
3. Abutment/ chronometric cap which is used. ROOTT R, ROOTT M/S/P- abutments, ROOTT C/CS/B/BS chronometric caps.
4. Prosthesis material (PEEK PMMA or Zirconium)

 Trate_ROOTT_External_platform_AntRot_peek_...	 Trate_ROOTT_External_platform_3mm
 Trate_ROOTT_External_platform_4mm	 Trate_ROOTT_External_platform_5mm
 Trate_ROOTT_External_platform_7mm	 Trate_ROOTT_External_platform_AntRot_Metal
 Trate_ROOTT_External_platform_Rot_peek_pm...	 Trate_ROOTT_External_platform_Rot_zro
 Trate_ROOTT_External_platform_TCE0_PCE0_pe...	 Trate_ROOTT_External_platform_TCE0_PCE0_zro
 Trate_ROOTT_External_platform_TCE1_PCE1_pe...	 Trate_ROOTT_External_platform_TCE1_PCE1_zro
 Trate_ROOTT_External_platform_TCE2_PCE2_pe...	 Trate_ROOTT_External_platform_TCE2_PCE2_zro
 Trate_ROOTT_External_platform_TCE3_PCE3_pe...	 Trate_ROOTT_External_platform_TCE3_PCE3_zro
 Trate_ROOTT_External_platform_TCES0_PCES0_...	 Trate_ROOTT_External_platform_TCES0_PCES0_...
 Trate_ROOTT_External_platform_TCES1_PCES1_...	 Trate_ROOTT_External_platform_TCES1_PCES1_...
 Trate_ROOTT_External_platform_TCES2_PCES2_...	 Trate_ROOTT_External_platform_TCES2_PCES2_...
 Trate_ROOTT_External_platform_TCEXS1_PCEXS...	 Trate_ROOTT_External_platform_TCEXS1_PCEXS...
 Trate_ROOTT_External_platform_TCEXS2_PCEXS...	 Trate_ROOTT_External_platform_TCEXS2_PCEXS...
 Trate_ROOTT_M_MP_Multi_Unit	 Trate_ROOTT_M_PCOM_peek_pmma
 Trate_ROOTT_M_PCOM_Rotational_peek_pmma	 Trate_ROOTT_M_PCOM_Rotational_zro
 Trate_ROOTT_M_PCOM_zro	 Trate_ROOTT_R_CRE_peek_pmma
 Trate_ROOTT_R_CRE_zro	 Trate_ROOTT_R_CRE07_peek_pmma
 Trate_ROOTT_R_CRE07_zro	 Trate_ROOTT_R_M1_Multi_Unit
 Trate_ROOTT_R_M1_PCOM_peek_pmma	 Trate_ROOTT_R_M1_PCOM_zro
 Trate_ROOTT_R_M1A15_PCOM_peek_pmma	 Trate_ROOTT_R_M1A45_PCOM_peek_pmma
 Trate_ROOTT_R_MS1_Multi_Unit	 Trate_ROOTT_R_MS1_PCOMS_peek_pmma
 Trate_ROOTT_R_MS1_PCOMS_zro	 Trate_ROOTT_R_PCO_peek_pmma
 Trate_ROOTT_R_PCO_zro	 Trate_ROOTT_R_PCO1_peek_pmma
 Trate_ROOTT_R_PCO1_zro	 Trate_ROOTT_R_PCO1S_peek_pmma
 Trate_ROOTT_R_PCO1S_zro	 Trate_ROOTT_R_PCO2_peek_pmma
 Trate_ROOTT_R_PCO2_zro	 Trate_ROOTT_R_PCO2S_peek_pmma
 Trate_ROOTT_R_PCO2S_zro	 Trate_ROOTT_R_PCO3_peek_pmma
 Trate_ROOTT_R_PCO3_zro	 Trate_ROOTT_R_PCO3S_peek_pmma
 Trate_ROOTT_R_PCO3S_zro	 Trate_ROOTT_R_PMAB
 Trate_ROOTT_S_Multi_Unit	 Trate_ROOTT_S_PCOMS_peek_pmma
 Trate_ROOTT_S_PCOMS_Rotational_peek_pmma	 Trate_ROOTT_S_PCOMS_Rotational_zro
 Trate_ROOTT_S_PCOMS_zro	

Instruction for Exocad using ROOTT R

Extraoral scan using PCO, PCOS

Step 1. Upload .slt file to Exocad software. Extraoral scan posts shall be visible.





Step 2. Choose abutment that will replace scan post.


Implant + abutment, material

e.g. ROOTT R + PCO1 zro

TRATE® ROOTT® - ROOTT R - PCO1 zro



**Intraoral or extraoral scan,
screw channel angle**
e.g. PCO1 25µm Extraoral
scan 0°



Cement gap size in microns, screw channel

e.g. PCO1 25µm screw
channel standard.

CRE.07

Short instruction on new version of CRE

It is important to know which version of CRE you are using, as they come with different measures. Mistakenly identified CRE parts can lead to inaccuracies and confusion, which can directly affect the work of implantologist and dental technicians.

The main disparity between CRE.07 and CRE.06 – height is of utmost importance during the processes of scanning while using CRE as a scanning body. Our digital library has individual options for each CRE and can be identified together with HE-END parts.

The main differences

CRE.07 is more streamlined, shorter and slightly narrower



CRE identification inside the blister

It is easy to identify the version of CRE when it is assembled with HE and the Implant:

CRE.06 height in this stand – 4,10mm

CRE.07 height in this stand – 2,76mm



CRE identification after insertion

CRE.06 height from the implant level is – 8,02mm


CRE.07 height from the implant level is – 6,68mm



Instruction for Exocad using CRE 06

Direct scan using CRE.06 USE CRE scan post with all digital ROOTT R digital ti base and abutments.


Step 1. Upload .slt file to Exocad software. Extraoral scan posts shall be visible.



Step 2. Choose abutment that will replace scan post.

Implant + abutment, material

e.g. ROOTT R + CRE zro



Abutment, Cement gap size in microns, screw channel (standard/ narrow)

e.g. PCO1 25µm screw channel standard.

Instruction for Exocad using CRE 07

Direct scan using CRE.07 USE CRE scan post with all digital ROOTT R digital ti base and abutments.

Step 1. Upload .slt file to Exocad software. Extraoral scan posts shall be visible.



Step 2. Choose abutment that will replace scan post.

Implant + abutment, material

e.g. ROOTT R + PCO1 zro

TRATE® ROOTT® - ROOTT R - PCO1 zro

CRE07 25µm Screw channel standard

CRE07 35µm Screw channel standard

CRE07 50µm Screw channel standard

CRE07 25µm Screw channel narrow

CRE07 35µm Screw channel narrow

CRE07 50µm Screw channel narrow

SELECT IMPLANT

Full library Favorites only

TRATE® ROOTT® - ROOTT R - ▼ ★

TRATE® ROOTT® - ROOTT R - CRE07 zro

TRATE® ROOTT® - ROOTT R - CRE peek pmma

TRATE® ROOTT® - ROOTT R - CRE zro

TRATE® ROOTT® - ROOTT R - M1 Multi - Unit

TRATE® ROOTT® - ROOTT R - M1 PCOM peek pmma

TRATE® ROOTT® - ROOTT R - M1 PCOM zro

TRATE® ROOTT® - ROOTT R - MS1 Multi - Unit

TRATE® ROOTT® - ROOTT R - MS1 PCOMS peek pmma

TRATE® ROOTT® - ROOTT R - MS1 PCOMS zro

Choose scanned detail e.g.
 CRE 25µm CRE 07 and
 screw channel angle 0°, 5°,
 10°, 15°, 20

CRE 25µm Extraoral scan 0°

CRE 25µm Intraoral scan 0°

CRE 25µm CRE 06 0°

CRE 25µm CRE 07 0°

CRE 25µm CRE 06 HE-END 0°

CRE 25µm CRE 07 HE-END 0°

CRE 25µm A1K 0°

Cement gap size in

microns, screw channel

e.g. PCO1 25µm screw
channel standard.

Step 3. Mark an area to detect Scan post position.

TRATE® ROOTT® - ROOTT R - MS1 PCOMS peek pmma

TRATE® ROOTT® - ROOTT R - PCOMS MS1 zro

TRATE® ROOTT® - ROOTT R - PCO1S peek pmma


TRATE® ROOTT® - ROOTT R - PCO1S zro

TRATE® ROOTT® - ROOTT R - PCO1 peek pmma


TRATE® ROOTT® - ROOTT R - PCO1 zro

TRATE® ROOTT® - ROOTT R - PCO2S peek pmma

TRATE® ROOTT® - ROOTT R - PCO2S zro

**ROOTT R**


Material - PEEK, PMMA or ZrO

**Necessary products to make a prosthesis**Abutment
PCOxScrews
SFPCOx/SFPCOxLScrewdriver
SDLB

Extraoral scan using PMAB

Step 1. Upload .stl file to Exocad software. Extraoral scan post SPCO shall be visible (p. 14, Step 1.)


Step 2. Choose Pre-milled abutment PMAB that will replace scan post SPCO.




Step 3. Create a framework and mill the abutment.



1.



2.



3.



4.

1. Individual abutment 2. Anatomical crown
3. PMAB Pre-milled abutment with analog AN 4. Anatomical framework and analog

Necessary products to make a prosthesis



Abutment
PMAB



Screws
S8/SL8



Screwdriver
SDLB



Crown+ custom
abutment

Intraoral scan

Step 1. Upload scanned model to Exocad software. Intraoral scan post shall be visible.



Step 2. Choose abutment that will replace scan post (see p. 11, Step 2.).

Implant + abutment, material

e.g. ROOTT R + MS1+ PCOMS zro

TRATE® ROOTT® - ROOTT R - MS1 PCOMS zro

SELECT IMPLANT

Full library Favorites only

TRATE® ROOTT® - ROOTT R - ▼ ★

PCOMS 35µm S▼ PCOMS 35µm Ir▼

Info about this library...

PCOMS 25µm Screw channel standard

PCOMS 35µm Screw channel standard

PCOMS 50µm Screw channel standard

PCOMS 25µm Screw channel narrow

PCOMS 35µm Screw channel narrow

PCOMS 50µm Screw channel narrow

Intraoral or extraoral scan, screw channel angle
e.g. PCOMS 35 µm SPCOMIOS MS1 scan 0°

PCOMS MS1 35µm SPCOMS MS1 scan 10°

PCOMS MS1 35µm SPCOMS MS1 scan 15°

PCOMS MS1 35µm SPCOMS MS1 scan 20°

PCOMS MS1 35µm SPCOMIOS MS1 scan 0°

PCOMS MS1 35µm SPCOMIOS MS1 scan 5°

PCOMS MS1 35µm SPCOMIOS MS1 scan 10°


PCOMS MS1 35µm SPCOMIOS MS1 scan 15°

PCOMS MS1 35µm SPCOMIOS MS1 scan 20°

Cement gap size in microns, screw channel size
e.g. PCOMS 35 µm Screw channel standard

ROOTT R Digital workflow

Digital workflow ROOTT R





ROOTT M ROOTT P

Scan posts



ROOTT M Scan posts correspond with multiunit platform connections. Scan posts indicate the exact position of the implant in the jaw. During the scanning process, the information about the position is transferred into digital format. According to Scan post position, height and direction, Scan-post is converted to abutment and sets the precise location of analogue. Get yourself acquainted with all types of ROOTT M Scan posts, which are developed for the precise manufacturing of bridge prosthesis.

	 SPCOM	 SPCOMIO	 SPCOMS	 SPCOMIOS
Extraoral scanning				
Intraoral scanning				
Reusable				
Multi-units				
Easily scannable				
Long body				
Short body				

Abutments

ROOTT M and ROOTT P digital library have two options for abutments. There are an abutment PCOM and a digital abutment MU that can be used directly with implant, the famous possibility to use Trate_ROOTT M P_Multi Unit (MU connection) is for metal frameworks .

ROOTT M ROOTT P abutment characteristics

- Angled access for tunnel from 0° to 20°;
- Only for bridges.
- Trate_ROOTT M P_Multi unit it's for direct connection
- Trate_ROOTT M_PCOM_ it's for cement or screw retained connections.

TRATE® ROOTT® - ROOTT M ▼ ★

TRATE® ROOTT® - ROOTT M P - Multi - Unit

TRATE® ROOTT® - ROOTT M - PCOM peek pmma

TRATE® ROOTT® - ROOTT M - PCOM zro

TRATE® ROOTT® - ROOTT R - CRE peek pmma

TRATE® ROOTT® - ROOTT R - CRE zro

TRATE® ROOTT® - ROOTT R - M1 PCOM peek pmma

TRATE® ROOTT® - ROOTT R - M1 PCOM peek pmma

TRATE® ROOTT® - ROOTT R - MS1 PCOMS peek pmma

TRATE® ROOTT® - ROOTT R - PCOMS MS1 zro



ROOTT M ROOTT P

Direct connection to implant, only in
Digital library

TRATE® ROOTT® - ROOTT M - PCOM zro

TRATE® ROOTT® - ROOTT R - CRE peek pmma

TRATE® ROOTT® - ROOTT R - CRE zro

TRATE® ROOTT® - ROOTT R - M1 Multi - Unit

TRATE® ROOTT® - ROOTT R - M1 PCOM peek pmma

TRATE® ROOTT® - ROOTT R - M1 PCOM zro

TRATE® ROOTT® - ROOTT R - MS1 Multi - Unit

TRATE® ROOTT® - ROOTT R - MS1 PCOMS peek pmma

TRATE® ROOTT® - ROOTT R - MS1 PCOMS zro




ROOTT M ROOTT P

PCOM abutmentMaterial - peek,
pmma or zro


PCOM/ PCOMS options

TRATE® ROOTT® - ROOTT M - PCOM Rotational peek pmma



Rotational connection geometry. Multi unit connection, between abutment and prosthesis.

TRATE® ROOTT® - ROOTT M - PCOM peek pmma





Connection geometry with rotation lock. Multi unit connection, between implant and abutment.

Instruction for Exocad using ROOTT M ROOTT P

Extraoral scan

Step 1. Scanned model upload to Exocad software. Extraoral scan posts shall be visible.



Step 2. Choose abutment that will replace scan post.

Implant + abutment, material

e.g. ROOTT M + PCOM zro

TRATE® ROOTT® - ROOTT M - PCOM zro

PCOM 25µm Screw channel standard

PCOM 35µm Screw channel standard

PCOM 50µm Screw channel standard

PCOM 25µm Screw channel narrow

PCOM 35µm Screw channel narrow

PCOM 50µm Screw channel narrow


SELECT IMPLANT

Full library Favorites only

TRATE® ROOTT® - ROOTT M ▼ ★

PCOM 25µm Scr▼ PCOM 25µm Ext▼

Info about this li



Click 'Best fit match'

Best fit n

Intraoral or extraoral scan, screw channel angle
e.g. PCOM 25 µm Extraoral scan 10°

PCOM 25µm Extraoral scan 0°

PCOM 25µm Extraoral scan 5°

PCOM 25µm Extraoral scan 10°

PCOM 25µm Intraoral scan 0°

PCOM 25µm Intraoral scan 5°


PCOM 25µm Intraoral scan 10°

PCOM 25µm Intraoral scan 15°


PCOM 25µm Intraoral scan 20°

PCOM 25µm Extraoral scan 15°


Cement gap size in microns, screw channel size
e.g. PCOM 25 µm screw channel standard
(Suitable for screw driver and screw)




Standard screw channels



Cement gap 25µm




Prosthesis in Exocad




Digital model with analog ANMD

Necessary products to make a prosthesis




Ti base
PCOM

+




Screw
SFPCOM

+




Screwdriver
SDLB

→




Created by Pesterev Evgeniy




Digital MU
abutment

+




Screw
SFPCOM

+



Screwdriver
SDLB

→



Created by Sokratis Agapiou



ROOTT S

Scan post



ROOTT S Scan-posts corresponds with a small multiunit platform connection. Scan-posts indicate the exact position of the implant in the jaw. During the scanning process, the information about the position is transferred into digital format. According to Scan-post position, height and direction, Scan-post is converted to abutment and sets the precise location of analogue. Get yourself acquainted with all types of ROOTT S Scan-posts, which are developed for the precise manufacturing of bridge prostheses.

*screw-retained restorations with a wide and secure fixation screw



SPCOMS scan post characteristic

- Extraoral scanning
- The long body allows comfortable usage and precise results of laboratory workflow
- Easily scannable
- Reusable



SPCOMIOS scan post characteristic


- Intraoral scanning
- The short body allows comfortable usage and precise results of implantologists' workflow
- Easily scannable
- Reusable

Abutments

ROOTT S digital library has two options for abutments. There are an abutment PCOMS and a digital abutment that can be used directly without a physical one.


ROOTT S abutments characteristics

- Angled access for tunnel from 0° to 20°;
- Only for bridges.



ROOTT S

Direct connection to implant, only in Digital library Material - peek, pmma or zro




ROOTT S

Material - peek, pmma or zro

Instruction for Exocad using ROOTT S

Extraoral scan

Step 1. Scanned model upload to Exocad software. Extraoral scan posts shall be visible.



Step 2. Choose abutment that will replace scan post.

Implant + abutment, material

e.g. ROOTT S + PCOMS peek pmma

TRATE® ROOTT® - ROOTT S - PCOMS peek pmma

Cement gap size in microns, screw channel
e.g. PCOMS 90 µm Screw channel standard

PCOMS 70µm Screw channel standard

PCOMS 90µm Screw channel standard

PCOMS 110µm Screw channel standard

PCOMS 70µm Screw channel narrow

PCOMS 90µm Screw channel narrow


PCOMS110µm Screw channel narrow

TRATE® ROOTT® - ROOTT S -▼

PCOMS 90µm S▼

PCOMS 90µm E▼

Info about this library...



Click 'Best fit matching' when done

Best fit matching

Intraoral or extraoral scan, screw channel angle.
e.g. PCOMS 90 µm Extraoral scan 0°

PCOMS 90µm Extraoral scan 0°

PCOMS 90µm Extraoral scan 5°

PCOMS 90µm Extraoral scan 10°

PCOMS 90µm Intraoral scan 0°



PCOMS 90µm Intraoral scan 5°

PCOMS 90µm Intraoral scan 10°


PCOMS 90µm Intraoral scan 15°

PCOMS 90µm Intraoral scan 20°

PCOMS 90µm Extraoral scan 15°

Step 4. Mark an area to detect Scan post position.**Intraoral scan****Step 1.** Upload scanned model to Exocad software. Intraoral scan post shall be visible.**Step 2.** Choose abutment that will replace scan post (see p. 7 , Step 2.).

Necessary products to make a prosthesis



What to use. All ROOTT M/S/P implants on six




Construction: Prosthesis on 8 ROOTT M,MS, MP implants and PCOM abutments



- TRATE® ROOTT® - ROOTT M P - Multi - Unit
- TRATE® ROOTT® - ROOTT S - PCOMS peek pmma
- TRATE® ROOTT® - ROOTT M - PCOM zro

Digital workflow **ROOTT M P S**



* MU abutment is only accessible in digital library with angulation option and used with SPCOMS screw for ROOTT S, SPCOM screw for ROOTT M/P.



ROOTT C ROOT CS ROOT B ROOT BS

Scanable superstructures



ROOTT C, CS, B, BS Scan-post give a wide range of options. There are 45 different ways to scan with intraoral and extraoral scanners – choose from TRA, HE TOEA, TOE, TOES, telescopic abutments or External platform varieties. Scan-posts have a few height options that open possibilities for different clinical cases and patients mouth.

Transfers can be used as scan-posts that make workflow more precise and effective. Scan-posts indicate the exact position and the depth of the screwed implant in the jaw. During the scanning process, the information about the position is transferred into digital format. According to the Scan-post position, Scan-post is converted to telescopic abutment and sets the precise location of analogue. Get yourself acquainted with one-piece abutments which are suitable even for complex clinical cases.



TRA



HE



TOEA



TOE



TOES

TCE, PCE
TCES, PCES
TCEXS, PCEXSExternal platform
7, 5, 4, 3 mm

Abutments

ROOTT C, CS, B, BS abutments are metal and plastic connectors between an implant and a crown. They do not require a screw – making it the best solution for time-saving and comfort. Both TCE and PCE are the same size; however, the material is different – titanium or peek. Plastic peek gives more amortization and cushioning in the tooth, making the bite more comfortable and reducing fracture risk.



TCE0
0 mm

TCE1
1 mm

TCE2
2 mm

TCE3
3 mm



TCEX0
0 mm

TCEX1
1 mm

TCEXS1
1 mm

TCEXS2
2 mm

Two ways to make direct prosthesis connection



ROOTT C Digital analog ANED, External platform




For Single crowns

TRATE® ROOTT® External platform - AntRot peek pmma

TRATE® ROOTT® External platform - AntRot_zro

In this way prosthesis milling with retention inside, which results- minimal prosthesis movement in milled zoned using an external platform.




For Bridges and arch

TRATE® ROOTT® External platform - Rot peek pmma

TRATE® ROOTT® External platform - Rot_zro

In this way, prosthesis milling is without any retention, resulting rotational prosthesis movement in milled zoned using an external platform.



Instruction for Exocad using ROOTT C CS B BS

Extraoral scan

Step 1. Scanned model upload to Exocad software. Extraoral scan posts shall be visible.



Step 2. Choose abutment that will replace scan post.

Abutment, material

e.g. External platform TCE1 PCE1 peek pmma

External platform - TCE1 PCE1 peek pmma

Full library Favorites only

TRATE® ROOTT® External pla...

TRATE® ROOTT® External platform - TCE1 PCE1 peek pmma

Scan post
e.g. TCE0 PCE0

TCE1 PCE1 90µm TOES Intraoral scan

TCE1 PCE1 90µm ANE Intraoral scan

TCE0 PCE0

TCE1 PCE1

TCE2 PCE2

TCE3 PCE3

TCEXS1 PCEXS1

TCEXS2 PCEXS2

Cement gap size in microns, screw channel
e.g. TCE PCE 70 µm screw channel standard
(Another possibility- no screw inside the crown)

TCE1 PCE1 70µm

TCE1 PCE1 90µm


TCE1 PCE1 110µm

TCE1 PCE1 70µm no screw channel


TCE1 PCE1 90µm no screw channel

TCE1 PCE1 110µm no screw channel


Step 3. Mark an area to detect Scan post position. Detected scan post shall change a color.



Step 4. Create a crown's anatomy, print a model and insert a digital analog in a digital model.




External platform with transparent caps.



External platform with cap's gap and





Digital analogs

Digital analogues indicate the exact position of the implant in the jaw. Therefore, if the scan post is not screwed properly, it could lead to an inaccurate position of digital analogue. Digital analogues could be used only with printed models.

ROOTT **R**



AND




AND analogs inserted in to the printed model

ROOTT **P**
ROOTT **M**



ANMD



ANMD Digital analogs suitable for ROOTT M/ ROOTT P

ROOTT **S**



ANMSD



ANED Digital analogs suitable for ROOTT C, ROOTT CS, ROOTT B, ROOTT BS

ROOTT **C**
ROOTT **CS**
ROOTT **B**
ROOTT **BS**



ANED




Digital analog AND inside the model



Analog shows the accurate implants position



Exocad let's change gap between model and digital analog




Analog ANED with empty space in model


Model Creator

An accurate printed model is needed to make a precise prosthesis. Therefore, it is crucial to know Model Creator software.

Note 1. Lift up a model in a way that analogue fits in it. If analogues are visible, the jaw model should be higher. If analogs are invisible, model is in the right position.



Incorrect- model is too low











































































Correct- good model position

PMAB milling opportunities

MILLING MACHINE			LINK TO HOLDER		
VHF	Available with library		Exocad	3Shape- Approved	
Roland DWX-42W	Available		Medentika Preface	Adapter	
ROLAND DWX-50	Available		Medentika Preface	Adapter	
IMES-ICORE CORITEC 250I	Available		loading		
Zfx	Available		loading		
LabTec	Available		Preface Adapter		
Dental direct S-Tec, K-tec	Available		loading		
Arum 400	Impossible		loading		

Libraries Folders sturcture


ROOTT R

	CRE		Analog AND		Scan posts SPCO SPCOIO		Screw SLT8		Connection geometry CRE.06_0.070_0bp CRE.06_0.110_0sp ...		Abutment CRE
ROOTT-R_TRATE_ROOTT-CRE											
	M1 Multi Unit		Analog AN		Scan posts SPCO SPCOIO SPCOM+M1 SPCOMIO+M1		Screw SFPCOM		Connection geometry PCOM.03_BRIDGE_0.025_15 PCOM.03_BRIDGE_0.025_20 ...		Digital analog MU
ROOTT-R_TRATE_ROOTT_M1-Multi-Unit											
	M1 PCOM		Analog AN		Scan posts SPCO SPCOIO SPCOM+M1 SPCOMIO+M1		Screw SFPCOM		Connection geometry PCOM.03_0.070_0bp PCOM.03_0.070_0sp PCOM.03_0.070_5bp ...		Abutment PCOM+M1
ROOTT-R_TRATE_ROOTT_M1-PCOM											
	MS1 Multi Unit		Analog AN		Scan posts SPCO SPCOIO SPCOMS+MS1 SPCOMIOS+MS1		Screw SFPCOMS		Connection geometry PCOMS.03_BRIDGE_0.025_15 PCOMS.03_BRIDGE_0.025_20 ...		Abutment PCOMS+MS1
ROOTT-R_TRATE_ROOTT_MS1-Multi-Unit											
	MS1 PCOMS		Analog AN		Scan posts SPCO SPCOIO SPCOMS+MS1 SPCOMIOS+MS1		Screw SFPCOMS		Connection geometry PCOMS.03_0.070_0bp ...		Abutment PCOMS+MS1
ROOTT-R_TRATE_ROOTT_MS1-PCOMS											
	PCO		Analog AN		Scan posts SPCO SPCOIO		Screw SFPCO1		Abutment PCO		Connection geometry PCO.04_0.070_0bp PCO.04_0.070_0sp PCO.04_0.070_5bp ...
ROOTT-R_TRATE_ROOTT_PCO											
	PCO1		Analog AN		Scan posts SPCO SPCOIO		Screw SFPCO1		Abutment PCO1		Connection geometry PCO1.04_0.070_0bp ...
ROOTT-R_TRATE_ROOTT_PCO1											
	PCO2		Analog AN		Scan posts SPCO SPCOIO		Screw SFPCO2		Abutment PCO2		Connection geometry PCO2.04_0.070_0bp ...
ROOTT-R_TRATE_ROOTT_PCO2											
	PCO3		Analog AN		Scan posts SPCO SPCOIO		Screw SFPCO3		Abutment PCO3		Connection geometry PCO3.04_0.070_0bp ...
ROOTT-R_TRATE_ROOTT_PCO3											
	PCO1S		Analog AN		Scan posts SPCO SPCOIO		Screw SFPCO1S		Abutment PCO1S		Connection geometry PCO1S.04_0.070_0bp ...
ROOTT-R_TRATE_ROOTT_PCO1S											
	PCO2S		Analog AN		Scan posts SPCO SPCOIO		Screw SFPCO2S		Abutment PCO2S		Connection geometry PCO2S.04_0.070_0bp ...
ROOTT-R_TRATE_ROOTT_PCO2S											
	PCO3S		Analog AN		Scan posts SPCO SPCOIO		Screw SFPCO3S		Abutment PCO3S		Connection geometry PCO3S.04_0.070_0bp ...
ROOTT-R_TRATE_ROOTT_PCO3S											
	PMAB		Analog AN		Scan posts SPCO SPCOIO		Screw SL8		Abutment PMAB		
ROOTT-R_TRATE_ROOTT_PMAB											












ROOTT M ROOTT P



ROOTT S



ROOTT C ROOTT CS ROOTT B ROOTT BS

		Analog ANE		Connection geometry TCE0.01_0.070 TCE0.01_0.090 ...	Abutment TCE0
	ROOTT-R_TRATE_ROOTT_TCE0-PCE0				
		Analog ANE		Connection geometry TCE1.01_0.070 TCE1.01_0.090 ...	Abutment TCE1
	ROOTT-R_TRATE_ROOTT_TCE1-PCE1				
		Analog ANE		Connection geometry TCE2.01_0.070 TCE2.01_0.090 ...	Abutment TCE2
	ROOTT-R_TRATE_ROOTT_TCE2-PCE2				
		Analog ANE		Connection geometry TCE3.01_0.070 TCE3.01_0.090 ...	Abutment TCE3
	ROOTT-R_TRATE_ROOTT_TCE3-PCE3				
		Analog ANE		Connection geometry TCES0.01_0.070 TCES0.01_0.090 ...	Abutment TCES0
	ROOTT-R_TRATE_ROOTT_TCES0-PCES0				
		Analog ANE		Connection geometry TCES1.01_0.070 TCES1.01_0.090 ...	Abutment TCES1
	ROOTT-R_TRATE_ROOTT_TCES1-PCES1				
		Analog ANE		Connection geometry TCES2.01_0.070 TCES2.01_0.090 ...	Abutment TCES2
	ROOTT-R_TRATE_ROOTT_TCES2-PCES2				
		Analog ANE		Connection geometry TCESX1.01_0.070 TCESX1.01_0.090 ...	Abutment TCESX1
	ROOTT-R_TRATE_ROOTT_TCESX1-PCEXS1				
		Analog ANE		Connection geometry TCESX2.01_0.070 TCESX2.01_0.090 ...	Abutment TCESX2
	ROOTT-R_TRATE_ROOTT_TCESX2-PCEXS2				