

# INSTRUCTION FOR IMPLANT REMOVAL

ROOT™ implant system has been thoroughly investigated, deeply clinically tested and has a big range of implant types, accessories and also, satisfied users.

Usually, a successful implantation depends on good surgical and prosthodontic planning. Especially using a 3D dental cone beam computed tomography planning and if it's possible producing an implant surgical guide. Also the presurgical patient preparation plays an important role.

It is known that errors can occur. Dental implantation failure can happen for different reasons:

- Peri-implantitis;
- Bone resorption;
- When it's impossible to have proper prosthesis (poor placement).

**Important: It is the physician's responsibility to evaluate the current level of osseointegration and apply the necessary system.**

Every implantologist has a situation in which immediately after the implantation doctor notices wrong positioning or angulation of the implant. Sometimes it can be more serious like perforation of the buccal or lingual bone cortex or dental implant in very close contact with the content of the mandibular canal, or dental implant that perforates sinus floor and etc. In this situation it is easy to remove the implant by using the same hand instrument used for implantation but in a counterclockwise direction. In situations like this, explantation is not a very hard operation, because there is only implant primary stability which is mechanical.

**Important: These recommended instruments can only be used for non-osseointegrated implants to be removed. It's the physician's responsibility to choose the right length of a removal tool.**

REF	Description
TW50	Torque wrench (10-50 Ncm)



DW	Direct wrench
----	---------------



## ROOTFORM

IT	Insertion tool, short
----	-----------------------



ITL	Insertion tool, long
-----	----------------------



## COMPRESSIVE, COMPRESSIVE S, BASAL, BASAL SS

ITES	Insertion tool, external platform, extra short
------	--



ITE	Insertion tool, external platform, short
-----	--



ITEL	Insertion tool, external platform, long
------	---



ITEXL	Insertion tool, external platform, extra long
-------	---



## COMPRESSIVE M, COMPRESSIVE MP

ITM0	Insertion tool for multiunit, height 0 mm
------	---



ITM10	Insertion tool for multiunit, height 10 mm
-------	--



ITM	Insertion tool for multiunit, short
-----	-------------------------------------



ITML	Insertion tool for multiunit, long
------	------------------------------------



ITMXL	Insertion tool for multiunit, extra long
-------	--



## COMPRESSIVE MS

ITMS	Insertion tool for small multiunit, short
------	---



ITMSL	Insertion tool for small multiunit, long
-------	--



ITMSXL	Insertion tool for small multiunit, extra long
--------	--

## COMPRESSIVE K

ITEKS	Insertion tool, K platform, extra short
-------	---



ITEK	Insertion tool, K platform, short
------	-----------------------------------




ITEKL	Insertion tool, K platform, long
-------	----------------------------------



ITEKXL	Insertion tool, K platform, extra long
--------	--



# Tools for implant removal




## Removal steps for non-osseointegrated two-piece and one-piece implants


- 1 Insert and lock the removal tool for a two-piece implant inside the implant, for one-piece implant - over implant.
- 2 Use a torque wrench TW50 and a direct wrench (DW) (for stabilisation). Then slowly unscrew the implant by applying force 50 Ncm (it is not recommended to use more than 100 Ncm force) to the removal tool and turning it counterclockwise.
- 3 Remove an implant with an insertion tool (pictures below).

**ATTENTION! Removed implants must not be reused!**

# Two-piece implant removal



# One-piece implant removal



# One-piece implant removal with multi-unit platform

