

# SURGICAL CASSETTE TRS-MINI

## instructions for use

The TRS-mini surgical cassette has a new design and layout:

- Compact and easy to carry
- Laser markings on the tray
- Has drill length check

For implant systems: **ROOTT R C C S M S**

How to open and close TRS-mini surgical cassette



Take a cassette so that the logo of ROOTT is facing towards you. Push the coverlid to the right side till every instrument in the cassette is visible and remove coverlid.




Put the coverlid with the ROOTT logo facing towards you in the grooves, and push it to the left side until it covers all of the instruments.

Surgical cassette must be cleaned, disinfected and sterilized before use and between each use. Do not use damaged instruments.

**Cassette does not include instruments, these are sold separately.**

# CONTENT



## Drills

**Lance drill DB2020** can be used for initial drilling by setting the drilling axis before using tapered drills. Drill has laser markings every 2 mm from 6 to 20 mm.



DB2020

**Drill D5506** is used for cortical drilling. As optional solution to widen the osteotomy.



D5506

**Tapered drills** are used for final cavity preparation. Use a drill with a diameter lower than the diameter of the implant to be placed. Drill to the appropriate depth, required for a specific case. Laser markings every 2 mm from 4 to 16 mm.



D2516

D3216

D3616

D4316

D5016

## Screwdriver


For all ROOTT dental implant system screws. Long screwdriver is suitable for all superstructures. Conical tip of the hex helps to grab screw.



SDL

## Implant drivers

**For external platform.** To control the implant insertion depth, have two lasered round lines at 2 and 4 mm.



Hole for lifting in case if tool is stacked

**For internal platform.** Have length indicators.



Length indicators 1 mm between each line

**For M platform via carrier.** Consists of two parts: implant driver and screwdriver.



ITHES

SITHES

## Abutment extractor

For easy superstructure removal in case if conical connection hold tight a part inside of ROOTT R implant.



SR

## Torque wrench

Torque wrench is suitable for all instruments with head for ratchet. Maximal torque is 50 Ncm.



TW50

## STEP 1: Preparing cavity


STEP 2: Implant insertion via carrier

STEP 3: Removing carrier

STEP 4: Implant insertion directly

### Example:

ROOTT R implant  
Ø 3.8 mm  
Length – 12 mm



- Step 1** Take a lance drill DB2020 to make the first mark on the bone.
- Step 2** Use drill length check to check drill length. Use laser marking to identify position with regards to intraosseous length of implant. Drill to the same depth as implant length.
- Step 3** Take a tapered drill D2516 to enlarge the diameter of the hole. Use wider tapered drills to prepare hole for wider diameter implants.

### \*Optional

If after using the previous drill the torque is still more than 50 Ncm while inserting the implant, the cavity has to be widened. Just widen the osteotomy with drill D5506.

### NOTE!

Not suitable for all implants and all bone densities. See table 1.


Full drilling  
protocol  
for ROOTT R  
implants



			
Implant	D4 BONE	D2-D3 BONE	D1 BONE
Ø 3.0 mm	X	DB2020 D2516	X
Ø 3.5 mm	DB2020 D2516	X	DB2020 D2516 D3216
Ø 3.8 mm	X	DB2020 D2516 D3216	DB2020 D2516 D3216 D3616
Ø 4.2 mm	X	DB2020 D2516 D3216	DB2020 D2516 D3216 D3616
Ø 4.8 mm	DB2020 D2516 D3216 D3616	X	DB2020 D2516 D3216 D3616 D4316
Ø 5.5 mm	X	DB2020 D2516 D3216 D3616 D4316 D5016	X

Table 1

**Example:**  
ROOTT R implant  
Ø 3.8 mm  
Length – 12 mm



**Step 1** Take an implant driver for external platform ITE/ITEL for inserting an implant via a carrier.



**Step 2** Insert ITE/ITEL into the torque wrench TW50 and tighten the implant by rotating the wrench clockwise.

When the set torque is reached, the scale sleeve snaps around the axis in the wrench head.  
The release can be heard and felt.



Do not continue to use the wrench after the torque is achieved.  
The wrench or dental components could be damaged.

**Example:**  
 ROOTT R implant  
 Ø 3.8 mm  
 Length – 12 mm



**Step 1** Take a multipurpose screwdriver SDL for screwing & unscrewing any screw of ROOTT dental implant system.

Due to the conical tip of the hex, it is more manageable to take out the screw from the superstructure. Therefore if struggling to remove the screw from the superstructure, movement side to side before pulling out is allowed.

**Step 2** Unscrew screw SLT8 from CRE.


**Step 3** Take an abutment extractor SR. Screw SR instead of your screw until part is released and remove CRE.

Abutment extractor SR - for easy superstructure removing in case if your conical connection hold tight a part inside of ROOTT R implant



## STEP 4: Implant insertion directly

**Example:**  
ROOTT R implant  
Ø 3.8 mm  
Length – 12 mm



**Step 1** Take an implant driver for internal platform IT/ITL for inserting ROOTT R implants when the carrier part is removed.

**Step 2** Place IT/ITL to torque wrench TW50 and insert implant to the prepared hole.

When the set torque is reached. The scale sleeve snaps around the axis in the wrench head. The release can be heard and felt.



Do not continue to use the wrench after the torque is achieved. The wrench or dental components could be damaged.

For detailed information  
please read ROOTT R  
implant placement  
protocol






**Example:**  
ROOTT C implant  
Ø 4.0  
Length – 12 mm



C4012



**Step 1** Take a lance drill DB2020 to make the first mark on the bone.

**Step 2** Use drill length check to check drill length. Use laser marking to identify position with regards to intraosseous length of implant. Drill to the same depth as implant length.

**Step 3** Take a tapered drill D2516 to enlarge the diameter of the hole. Use wider tapered drills to prepare hole for wider diameter implants.

*\*Optional*

If after using the previous drill the torque is still more than 50 Ncm while inserting the implant, the cavity has to be widened. Just widen the osteotomy with drill D5506.

### NOTE!

Not suitable for all implants and all bone densities. See table 2.


Full drilling  
protocol for  
ROOTT C/CS  
implants



 Implant	D4 BONE		D2-D3 BONE		D1 BONE	* - drill to the depth as specified  *C3006, C3008 - 4mm C3010, C3012 - 6 mm C3014, C3016 - 8 mm C3016, C3020 -10 mm  *C3506, C3508 - 4mm C3510, C3512 - 6 mm C3514, C3516 - 8 mm C3516, C3520 -10 mm  *C4006, C4008 - 4mm C4010, C4012 - 6 mm C4014, C4016 - 8 mm C4016, C4020 -10 mm  *C4506, C4508 - 4mm C4510, C4512 - 6 mm C4514, C4516 - 8 mm C4516, C4520 -10 mm  *C5006, C5008 - 4mm C5010, C5012 - 6 mm C5014 - 8 mm  *C5506, C5508 - 4mm C5510, C5512 - 6 mm C5514 - 8 mm
	DB2020		X		DB2020 D2516*	
	X		DB2020 D2516*		X	
	6-20 mm DB2020 D2516*		X		6-8 mm DB2020 D2516 D3216*	
	X		6-20 mm DB2020 D2516 D3216*		6-20 mm DB2020 D2516 D3216* D3616*	
	6-8 mm DB2020 D2516 D3216	10-14 mm DB2020 D2516* D3216*	6-8 mm DB2020 D2516 D3216 D3616	10-14 mm DB2020 D2516 D3216* D3616*	X	
	X		6-8 mm DB2020 D2516 D3216 D3616 D4316*	10-14 mm DB2020 D2516 D3216* D3616* D4316*	X	

Table 2

**Example:**  
ROOTT C implant  
Ø 4.0  
Length – 12 mm



**Step 1** Take an implant driver for external platform ITE/ITEL for inserting an implant via a carrier.

**Step 2** Place ITE/ITEL to torque wrench TW50 and insert implant to the prepared hole.

When the set torque is reached, the scale sleeve snaps around the axis in the wrench head. The release can be heard and felt.



Do not continue to use the wrench after the torque is achieved. The wrench or dental components could be damaged.

For detailed information  
please read ROOTT C/CS  
implant placement  
protocol





#### Example:

ROOTT M implant

Ø 4.0

Length – 12 mm



**Step 1** Take a lance drill DB2020 to make the first mark on the bone.

**Step 2** Use drill length check to check drill length. Use laser marking to identify position with regards to intraosseous length of implant. Drill to the same depth as implant length.

**Step 3** Take a tapered drill D2516 to enlarge the diameter of the hole. Use wider tapered drills to prepare hole for wider diameter implants.

#### \*Optional

If after using the previous drill the torque is still more than 50 Ncm while inserting the implant, the cavity has to be widened. Just widen the osteotomy with drill D5506.

#### NOTE!

Not suitable for all implants and all bone densities. See table 3.

Full drilling  
protocol for  
ROOTT M  
implants





Implant	D4 BONE	D2-D3 BONE	D1 BONE	* - drill to the depth as specified
Ø 3.0 mm	DB2020	X	DB2020 D2516*	*C3008m - 4mm C3010m, C3012m - 6 mm C3014m, C3016m - 8 mm C3016m, C3020m -10 mm
Ø 3.5 mm	X	DB2020 D2516*	X	*C3506m, C3508 m- 4mm C3510m, C3512m - 6 mm C3514m, C3516m - 8 mm C3516m, C3520m -10 mm
Ø 4.0 mm	X	DB2020 D2516* D3216*	DB2020 D2516 D3216* D3616*	*C4006, C4008 - 4mm C4010, C34012 - 6 mm C4014, C4016 - 8 mm
Ø 5.0 mm	X	DB2020 D2516 D3216* D3616* D4316*	X	*C5006, C5008 - 4mm C5010, C5012 - 6 mm C5014 - 8 mm
Ø 6.0 mm	X	DB2020 D2516 D3216 D3616* D4316* D5016*	X	*C6006, C6008 - 4mm C6010, C6012 - 6 mm C6014 - 8 mm

Table 3


**Example:**  
 ROOTT M implant  
 Ø 4.0  
 Length – 12 mm



C4012m



ITHES  
 + SITHES



**Step 1** Take an implant driver ITHES for inserting ROOTT M implants via the carrier.

**Step 2** Place ITHES to torque wrench TW50 and insert implant to the prepared hole.



When the set torque is reached, the scale sleeve snaps around the axis in the wrench head. The release can be heard and felt.



Do not continue to use the wrench after the torque is achieved. The wrench or dental components could be damaged.

**Example:**

ROOTT M implant  
Ø 4.0  
Length – 12 mm



**Step 1** Carrier can be removed without removing the implant driver ITHES using its screwdriver SITHES.

Use screwdriver SITHES to unscrew SHFCM screw and remove carrier HFCM.

**Step 2** Take a multipurpose screwdriver SDL for screwing & unscrewing any screw of ROOTT dental implant system.

Due to the conical tip of the hex, it is more manageable to take out the screw from the superstructure. Therefore if struggling to remove the screw from the superstructure, movement side to side before pulling out is allowed.

**Step 3** Unscrew screw SHFCM and remove carrier HFCM.

For detailed information  
please read ROOTT M  
implant placement  
protocol





#### Example:

ROOTT S implant

Ø 3.0

Length – 12 mm



C3012ms



DB2020




D2516



D5506

\*Optional



**Step 1** Take a lance drill DB2020 to make the first mark on the bone.

**Step 2** Use drill length check to check drill length. Use laser marking to identify position with regards to intraosseous length of implant. Drill to the same depth as implant length.

**Step 3** Take a tapered drill D2516 to enlarge the diameter of the hole. Use wider tapered drills to prepare hole for wider diameter implants.

*\*Optional*

If after using the previous drill the torque is still more than 50 Ncm while inserting the implant, the cavity has to be widened. Just widen the osteotomy with drill D5506.

#### NOTE!

Not suitable for all implants and all bone densities. See table 4.

Full drilling  
protocol for  
ROOTT S  
implants





Implant	D4 BONE	D2-D3 BONE	D1 BONE	* - drill to the depth as specified
Ø 3.0 mm	DB2020	X	DB2020 D2516*	*C3008ms - 4mm C3010ms, C3012ms - 6 mm C3014ms, C3016ms - 8 mm
Ø 3.5 mm	X	DB2020 D2516*	X	*C3506ms, C3508ms- 4mm C3510ms, C3512ms - 6 mm C3514ms, C3516ms - 8 mm

Table 4

**Example:**  
 ROOTT S implant  
 Ø 3.0  
 Length – 12 mm




C3012ms



ITE

ITEL



**Step 1** Take an implant driver for external platform ITE/ITEL for inserting an implant via the carrier.

**Step 2** Place ITE/ITEL to torque wrench TW50 and insert implant to the prepared hole.



When the set torque is reached, the scale sleeve snaps around the axis in the wrench head. The release can be heard and felt.



Do not continue to use the wrench after the torque is achieved. The wrench or dental components could be damaged.

**Example:**

ROOTT S implant  
Ø 3.0  
Length – 12 mm



**Step 1** Take a multipurpose screwdriver SDL for screwing & unscrewing any screw of ROOTT dental implant system.

Due to the conical tip of the hex, it is more manageable to take out the screw from the superstructure. Therefore if struggling to remove the screw from the superstructure, movement side to side before pulling out is allowed.

**Step 2** Unscrew screw SCREM and remove carrier CREM.

For detailed information  
please read ROOTT M  
implant placement  
protocol





**ROOTT**

TRATE AG  
Bahnhofstrasse 16  
6037 Root  
Switzerland

+41 41 450 01 01  
info@trate.com  
trate.com

**TRATE**

Medical devices under these instructions are in compliance with established in EU regulatory requirements.



TRATE UAB  
Jonavos st. 254  
44110 Kaunas  
Lithuania



TRATE AG  
Bahnhofstrasse 16  
6037 Root  
Switzerland