


SURGICAL CASSETTE TRS

instructions for use



- Easy opening cover
- Laser markings on the tray
- Stable-fitted tools for safe transportation
- Convenient bundle of tools
- Color-coded for better recognition
- Drills arranged from left to right in order of increasing diameter
- Drill length check
- The insert is easily removable


For implant systems: **ROOTT R C C S M**



TRS cassette's insert is designed to hold drills, implant drivers, screwdrivers and some other instruments. Can be easily removed.

Cassette does not include instruments, these are sold separately.

CONTENT




Drills

Twist drill D2016 can be used as pilot drills. These drills have a stopper, which stops the drill from making deeper hole than needed. Depth marks for easy visual reference while drilling or x-ray control. Have laser markings from 4 mm each 2 mm.



Tapered drills are used for cavity preparation. Use a drill with a diameter lower than the diameter of the implant to be placed. Drill to the appropriate depth, required for a specific case. Laser markings every 2 mm from 4 to 16 mm.



Handle

For supporting implant drivers, that are used with a wrench



Torque wrench

Torque wrench is suitable for all instruments with head for ratchet. Maximal torque is 70 Ncm.



Drill stops

The drill stop restricts drilling deeper than the predefined depth, to be used in sensitive indications to avoid the mandibular nerve or sinus floor. All drill stops can be mounted and dismantled manually. Color-coding allows easy identification of diameter size.



Implant drivers

For external platform. To control the implant insertion depth, have two lasered round lines at 2 and 4 mm. For ratchet and handpiece.



37 ITE

Hole for lifting in case if tool is stacked



38 ITEL

0²4



39 ITEH



For internal platform. Have length indicators.



40 IT

Length indicators 1 mm between each line



41 ITL



42 ITHL



For M platform (Viper). Consists of two parts: implant driver and screwdriver.



43 ITM10



44 ITM

For handpiece, coming soon

45 ITHEH



Screwdrivers

For all ROOTT dental implant system screws. Long screwdriver is suitable for all superstructures. Conical tip of the hex helps to grab screw.



46 SD



47 SDL




Abutment extractor

For easy superstructure removal in case if conical connection hold tight a part inside of ROOTT R implant.



48




Optional places

49

For additional drills, handpiece instruments, ratchet instruments.

Example:
 ROOTT R implant
 Ø 3.8 mm
 Length – 10 mm




- Step 1** Take a twist pilot drill D2016 to define the direction of the implant and to enlarge the diameter of the hole.
 *For long implants Rxx18 and Rxx20 should be used longer twist drills D2018 or .
- Step 2** Use drill length check to check drill length. Use laser marking to identify position with regards to intraosseous length of implant. Drill deeper 0.5 mm then the implant length for subcrestal position.
- Step 3** Take a tapered drill D2516 to enlarge the diameter of the hole. Use wider tapered drills to prepare hole for wider diameter implants.
- Optional step**** Use drill stops to control drilling depth.

Drilling protocol for ROOTT R implants or see Table 1

Implant	D4 BONE	D2-D3 BONE	D1 BONE
Ø 3.0 mm	D2016	D2016 D2516	D2016 D2516 D2816
Ø 3.5 mm	D2016 D2516	D2016 D2516 D2816	D2016 D2516 D2816 D3216
Ø 3.8 mm	D2016 D2516 D2816	D2016 D2516 D2816 D3216	D2016 D2516 D2816 D3216 D3616
Ø 4.2 mm	D2016 D2516 D2816 D3216	D2016 D2516 D2816 D3216 D3616	D2016 D2516 D2816 D3216 D3616 D4016
Ø 4.8 mm	D2016 D2516 D2816 D3216 D3616 D4016	D2016 D2516 D2816 D3216 D3616 D4016 D4316	D2016 D2516 D2816 D3216 D3616 D4016 D4316 D4616
Ø 5.5 mm	D2016 D2516 D2816 D3216 D3616 D4016 D4316 D4616	D2016 D2516 D2816 D3216 D3616 D4016 D4316 D4616 D5016	D2016 D2516 D2816 D3216 D3616 D4016 D4316 D4616 D5016 D5316



Example:
 ROOTT R implant
 Ø 3.8 mm
 Length – 10 mm



- Step 1** Take an implant driver for internal platform IT/ITL/ITHL and insert to ROOTT R implant.
- Step 2** Twist off breaking point – hold the instrument tight with one hand and with the other hold the holder, twisting until the implant separates from the holder.
- Step 3** Use implant driver IT/ITL with torque wrench TW70 to insert implant to prepared hole. Use implant driver ITHL with contra–angle handpiece to insert implant to prepared hole.
- Step 4** Use handle for implant driver DW for more precise implant insertion and to avoid glove grab.


Implant insertion via direct insertion – never exceed 100 Ncm.

Example:
 ROOTT R implant
 Ø 3.8 mm
 Length – 10 mm



- Step 1** Use a multipurpose screwdriver SD or SDL for screwing & unscrewing any screw of ROOTT dental implant system.
 Due to the conical tip of the hex, it is more manageable to take out the screw from the superstructure. Therefore if struggling to remove the screw from the abutment, movement side to side before pulling out is allowed.
 Never exceed the recommended tightening torque – 15 Ncm for the screw. Overtightening of abutment may lead to a screw fracture.
- Step 2** Use an abutment extractor SR for easier removing ROOTT R healing abutment or D2020 abutment. Screw SR instead of your screw until part will not be released and remove CRE. Abutment extractor SR – for easy superstructure removal in case if your conical connection holds tight a part inside of ROOTT R implant.

Example:
 ROOTT C implant
 Ø 4.0 mm
 Length – 10 mm




- Step 1** Take a twist pilot drill D2016 to define the direction of the implant and to enlarge the diameter of the hole.
 *For long implants Cxx18 and Cxx20 should be used longer twist drills D2018 or D2020.
- Step 2** Use drill length check to check drill length. Use laser marking to identify position with regards to intraosseous length of implant. Drill to the required length (see Table 2 or drilling protocol for ROOTT C/CS).
- Step 3** Take a tapered drill D2516 to enlarge the diameter of the hole. Use wider tapered drills to prepare hole for wider diameter implants.
- Optional step**** Use drill stops to control drilling depth.

Drilling protocol for ROOTT C/CS implants see Table 2

Implant	D4 bone		D2-D3 bone		D1 bone		*drill to depth as specified
Ø 3.0 mm	D2016 ↻		D2016*		D2016 D2516*		*C3006, C3008 – 4 mm C3010, C3012 – 6 mm C3014, C3016 – 8 mm C3018, C3020 – 10 mm
Ø 3.5 mm	D2016		D2016 D2516*		D2016 D2516* D2816*		*C3506, C3508 – 4 mm C3510, C3512 – 6 mm C3514, C3516 – 8 mm C3518, C3520 – 10 mm
Ø 4.0 mm	6-20 mm D2016* D2516*		6-20 mm D2016 D2516* D2816*		6-8 mm D2016 D2516 D2816* D3216*	10-20 mm D2016 D2516* D2816* D3216*	*C4006, C4008 – 4 mm C4010, C4012 – 6 mm C4014, C4016 – 8mm C4018, C4020 – 10 mm
Ø 4.5 mm	6-20 mm D2016 D2516* D2816*		6-20 mm D2016 D2516 D2816* D3216*		6-8 mm D2016 D2516 D2816 D3216* D3616*	10-20 mm D2016 D2516 D2816* D3216* D3616*	*C4506, C4508 – 4 mm C4510, C4512 – 6 mm C4514, C4516 – 8 mm C4518, C4520 – 10 mm
Ø 5.0 mm	6-8 mm D2016 D2516 D2816 D3216	10-14 mm D2016 D2516* D2816* D3216*	6-8 mm D2016 D2516 D2816 D3216 D3616	10-14 mm D2016 D2516 D2816* D3216* D3616*	6-8 mm D2016 D2516 D2816 D3216 D3616 D4016	10-14 mm D2016 D2516 D2816 D3216* D3616* D4016*	*C5006, C5008 – 4 mm C5010, C5012 – 6 mm C5014 – 8 mm
Ø 5.5 mm	6-8 mm D2016 D2516 D2816 D3216 D3616 D4016*	10-14 mm D2016 D2516* D2816* D3216* D3616* D4016*	6-8 mm D2016 D2516 D2816 D3216 D3616 D4016 D4316*	10-14 mm D2016 D2516 D2816* D3216* D3616* D4016* D4316*	6-8 mm D2016 D2516 D2816 D3216 D3616 D4016 D4316 D4616*	10-14 mm D2016 D2516 D2816 D3216* D3616* D4016* D4316* D4616*	*C5506, C5508 – 4 mm C5510, C5512 – 6 mm C5514 – 8 mm
Ø 6.5 mm	6-8 mm D2016 D2516 D2816 D3216 D3616 D4016 D4316 D4616*	10-14 mm D2016 D2516 D2816 D3216 D3616 D4016 D4316* D4616*	6-8 mm D2016 D2516 D2816 D3216 D3616 D4016 D4316 D4616 D5016	10-14 mm D2016 D2516 D2816 D3216 D3616 D4016 D4316 D4616* D5016*	6-8 mm D2016 D2516 D2816 D3216 D3616 D4016 D4316 D4616 D5016 D5316	10-14 mm D2016 D2516 D2816 D3216 D3616 D4016 D4316 D4616 D5016* D5316*	*C6506, C6508 – 4 mm C6510, C6512 – 6 mm


Example:
 ROOTT C implant
 Ø 4.0 mm
 Length – 10 mm



- Step 1** Take an implant driver for external platform ITE/ITEL/ITEH and put on ROOTT C/CS implant.
- Step 2** Twist off breaking point – hold the instrument tight with one hand and with the other hold the holder, twisting until the implant separates from the holder.
- Step 3** Use implant driver ITE/ITEL with torque wrench TW70 to insert implant to prepared hole. Use implant driver ITEH with a contra-angle handpiece to insert the implant to the prepared hole.
- Step 4** Use handle for implant driver DW for more precise implant insertion and to avoid glove grab.

Implant insertion via direct insertion:
 – from diameter 3.0 mm never exceed 117 Ncm
 – from diameter 3.5 mm never exceed 133 Ncm
 – from diameter 4.5-5.0 mm never exceed 238 Ncm
 – from diameter 5.5 mm never exceed 298 Ncm

Example:
 ROOTT M implant
 Ø 4.0 mm
 Length – 10 mm




- Step 1** Take a twist pilot drill D2016 to define the direction of the implant and to enlarge the diameter of the hole.
 *For long implants Cxx18m and Cxx20m should be used longer twist drills D2018 or D2020.
- Step 2** Use drill length check to check drill length. Use laser marking to identify position with regards to intraosseous length of implant. Drill to the required length (see table 3 or Drilling protocol for ROOTT M).
- Step 3** Take a tapered drill D2516 to enlarge the diameter of the hole. Use wider tapered drills to prepare hole for wider diameter implants.
- Optional step**** Use drill stops to control drilling depth.

Drilling protocol for ROOTT M implants see Table 3

Implant	D4 BONE		D2-D3 BONE		D1 BONE
Ø 2.5 mm	D2016 ↻		D2016		D2016
Ø 3.0 mm	D2016		D2016 D2516		D2016 D2516 D2816
Ø 3.5 mm	6-8 mm D2016 D2516	10-14 mm D2016 D2516 D2816	6-8 mm D2016 D2516 D2816	10-14 mm D2016 D2516 D2816 D3216	D2016 D2516 D2816 D3216
Ø 4.0 mm	6-8 mm D2016 D2516 D2816	10-14 mm D2016 D2516 D2816 D3216	6-8 mm D2016 D2516 D2816 D3216	10-14 mm D2016 D2516 D2816 D3216 D3616	D2016 D2516 D2816 D3216 D3616
Ø 5.0 mm	6-8 mm D2016 D2516 D2816 D3216 D3616 D4016	10-14 mm D2016 D2516 D2816 D3216 D3616 D4016 D4316	6-8 mm D2016 D2516 D2816 D3216 D3616 D4016 D4316	10-14 mm D2016 D2516 D2816 D3216 D3616 D4016 D4316 D4616	D2016 D2516 D2816 D3216 D3616 D4016 D4316 D4616


Example:
 ROOTT M implant
 Ø 4.0 mm
 Length – 10 mm



- Step 1** Take an implant driver for multi unit platform ITM10/ITM and put on ROOTT M implant. Carefully screw to the implant.
- Step 2** Twist off breaking point – hold the instrument tight with one hand and with the other hold the holder, twisting until the implant separates from the holder.
- Step 3** Use implant driver ITM10/ITM with torque wrench TW70 to insert implant to prepared hole.
- Step 4** Use handle for implant driver DW for more precise implant insertion and to avoid glove grab.

Implant insertion via direct insertion with insertion tools never exceed 133 Ncm.

Example:
 ROOTT M implant
 Ø 4.0 mm
 Length – 10 mm



Step 1 Use a multipurpose screwdriver SD or SDL for screwing & unscrewing any screw of ROOTT dental implant system.
 Due to the conical tip of the hex, it is more manageable to take out the screw from the superstructure. Therefore if struggling to remove the screw from the abutment, movement side to side before pulling out is allowed.

Never exceed the recommended tightening torque – 15 Ncm for the screw. Overtightening of abutment may lead to a screw fracture.



ROOTT

TRATE AG
Bahnhofstrasse 16
6037 Root
Switzerland

+41 41 450 01 01
info@trate.com
trate.com

TRATE

Medical devices under these instructions are in compliance with established in EU regulatory requirements.

 TRATE AG
Bahnhofstrasse 16
6037 Root
Switzerland

 TRATE UAB
Kauno m. sav. Kauno m.
Jonavos g 254
Kaunas 44110 Lithuania