



ROOTT

Creating restoration without screw channel
on External platform

guide for Exocad software



Step 1.

Close DentalCADApp if it is open.

Step 2.

Go to: DentalCADApp/config folder on your computer.

Step 3.

Make a backup copy of settings.xml file.

Step 4.

Open settings.xml file and enter the tag line:
<ImplantCloseHoles>true</ImplantCloseHoles>
between tag lines:
<ApplicationSettingsProviderCAD>
and
</ApplicationSettingsProviderCAD>

Step 5.

Save and close settings.xml file.

Note:

To turn screw channel hole feature back on, simply delete added tag line or change the value “true” to “false”.