



implant digital library guide for
Exocad software

Step 1.

Select the abutment of the implant.

The screenshot shows a software interface titled "Detect Implant Position". Below the title is a section labeled "Select library parts" which contains a dropdown menu. The dropdown is open, displaying a list of implant abutment options. The options are: Trate® ROOTT® C, CS, B, BS External platform height: 3mm; Trate® ROOTT® C, CS, B, BS External platform height: 4mm; Trate® ROOTT® C, CS, B, BS External platform height: 5mm; Trate® ROOTT® C, CS, B, BS External platform height: 7mm; Trate® ROOTT® M, P - PCOM; Trate® ROOTT® R - PCO; Trate® ROOTT® R - PCO1; Trate® ROOTT® R - PCOS; Trate® ROOTT® R - PCOR; and Trate® ROOTT® S - Multi-Unit. At the bottom of the dialog are two buttons: "BACK" and "NEXT".

Detect Implant Position

Select library parts

Trate® ROOTT® C, CS, B, BS External platform height: 3mm

Trate® ROOTT® C, CS, B, BS External platform height: 4mm

Trate® ROOTT® C, CS, B, BS External platform height: 5mm

Trate® ROOTT® C, CS, B, BS External platform height: 7mm

Trate® ROOTT® M, P - PCOM

Trate® ROOTT® R - PCO

Trate® ROOTT® R - PCO1

Trate® ROOTT® R - PCOS

Trate® ROOTT® R - PCOR

Trate® ROOTT® S - Multi-Unit

BACK NEXT

Step 2.

Select the cement gap, screw channel type and production material.

The screenshot shows the same "Detect Implant Position" dialog box. The "Select library parts" dropdown menu is now closed, and the selected option is "Trate® ROOTT® R - PCO1S". Below the dropdown is another dropdown menu labeled "Select an option". This second dropdown is open, showing a list of options for cement gap, screw channel type, and production material. The options are: 25µm Screw channel standard zro; 35µm Screw channel standard zro; 50µm Screw channel standard zro; 25µm Screw channel narrow zro; 35µm Screw channel narrow zro; 50µm Screw channel narrow zro; 70µm Screw channel standard peek/pmma; 90µm Screw channel standard peek/pmma; and 110µm Screw channel standard peek/pmma. The "BACK" and "NEXT" buttons are still visible at the bottom.

Detect Implant Position

Select library parts

Trate® ROOTT® R - PCO1S

Select an option

25µm Screw channel standard zro

35µm Screw channel standard zro

50µm Screw channel standard zro

25µm Screw channel narrow zro

35µm Screw channel narrow zro

50µm Screw channel narrow zro

70µm Screw channel standard peek/pmma


90µm Screw channel standard peek/pmma

110µm Screw channel standard peek/pmma

BACK NEXT

Step 3.

Select the scan post used for scan and screw channel angle.



Step 4.

Detect implant position by clicking marked point on scan and use “Best fit matching” function.

