



# ROOTT

Creating restoration without screw channel  
on External platform  
guide for Exocad software



**Step 1.**

Close DentalCADApp if it is open.

**Step 2.**

Go to: DentalCADApp/config folder on your computer.

**Step 3.**

Make a backup copy of settings.xml file.

**Step 4.**

Open settings.xml file and enter the tag line:  
<ImplantCloseHoles>true</ImplantCloseHoles>  
between tag lines:  
<ApplicationSettingsProviderCAD>  
and  
</ApplicationSettingsProviderCAD>

**Step 5.**

Save and close settings.xml file.

**Tip:**

To turn screw channel hole feature back on, simply delete added tag line or change the value “true” to “false”.