

# Digital library compatibility table



\* Material defines the range of cement gaps  
 \*\* Used only with ROOTT M  
 \*\*\* Possible abutment rotational change



Abutment	Screw	Material*	Gap between abutment and the crown (µm)	Library selection names	AMS	SPCOMS	SPCOMIOS	CREM+ HE_END	AM	SPCOM	SPCOMIO	PHFCM**	Direct Scan	Rotational option***	Screw channel options			
															Angulated 0-20°	Straight 0°	Standard	Narrow
	Sold separately (SFPSCOM)	Metal	0,15,20,25,35	Trate® ROOTT® S - Multi - Unit	●	●	●	●	—	—	—	—	●	Fully	●	●	●	—
	PCOMS PCOMS	Zirconium (zro) PEEK/PMMA	Zirconium 25, 35, 50 PEEK/PMMA 70, 90, 110	Trate® ROOTT® S - PCOMS Rotational Trate® ROOTT® S - PCOMS	●	●	●	●	—	—	—	—	●	Fully	—	●	●	●
	Sold separately (SFPSCOM)	Metal	0,15,20,25,35	Trate® ROOTT® M, P - Multi - Unit	—	—	—	—	●	●	●	●	●	Fully	●	●	●	—
	PCOM PCOM	Zirconium (zro) PEEK/PMMA	Zirconium 25, 35, 50 PEEK/PMMA 70, 90, 110	Trate® ROOTT® M, P - PCOM Rotational Trate® ROOTT® M, P - PCOM	—	—	—	—	●	●	●	●	●	Fully	—	●	●	●