




*** Possible abutment rotational change

																			Screw channel options			
																			Angulated 0-20°	Straight 0°	Standard	Narrow
	Abutment	Screw	Material*	Gap between abutment and the crown (µm)	Library selection names	Extraoral (SPCO)	Intraoral (SPCIO.02)	CRE.06	+HE_END	CRE.07	+HE_END	SPCOM +Mx(Ax)	SPCOMIO +Mx(Ax)	SPCOMS +MSx	SPCOMIOS +MSx	A1K	Rotational option***					
	M1+PCOM	Included	Zirconium (zro)	Zirconium 25, 35, 50	Trate® ROOTT® R - M1+PCOM	●	●	●	●	●	●	●	●	●	●	●	Fully rotational	●	●	●	●	
	MS1+PCOMS		Trate® ROOTT® R - MS1+PCOMS	●	●	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	
	M2+PCOM		PEEK/PMMA	PEEK/PMMA 70, 90, 110	Trate® ROOTT® R - M2+PCOM	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	●
	MS2+PCOMS			Trate® ROOTT® R - MS2+PCOMS	●	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	
	M3+PCOM			Trate® ROOTT® R - M3+PCOM	●	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	●
	MS3+PCOMS			Trate® ROOTT® R - MS3+PCOMS	●	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	●
	M4+PCOM			Trate® ROOTT® R - M4+PCOM	●	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	●
	MS4+PCOMS			Trate® ROOTT® R - MS4+PCOMS	●	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	●
	M1 multi-unit	Sold separetly (SFPCOM)	Metal	0, 15, 25, 35	Trate® ROOTT® R - M1 Multi - Unit	●	●	●	●	●	●	●	●	●	●	●		●	●	●	—	
	MS1 multi-unit				Trate® ROOTT® R - MS1 Multi - Unit	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	—
	M2 multi-unit				Trate® ROOTT® R - M2 Multi - Unit	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	—
	MS2 multi-unit				Trate® ROOTT® R - MS2 Multi - Unit	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	—
	M3 multi-unit				Trate® ROOTT® R - M3 Multi - Unit	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	—
	MS3 multi-unit				Trate® ROOTT® R - MS3 Multi - Unit	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	—
	M4 multi-unit				Trate® ROOTT® R - M4 Multi - Unit	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	—
	MS4 multi-unit				Trate® ROOTT® R - MS4 Multi - Unit	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	—
	Multi-unit angled	Included	Zirconium (zro)	Zirconium 25, 35, 50	TRATE® ROOTT® R - M1A15+PCOM	●	●	●	●	●	●	●	●	●	●	●	No	●	●	●	●	
	M1A15+PCOM		TRATE® ROOTT® R - M1A30+PCOM	●	●	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	
	M1A30+PCOM		TRATE® ROOTT® R - M1A45+PCOM	●	●	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	
	M1A45+PCOM		TRATE® ROOTT® R - M2A15+PCOM	●	●	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	
	M2A15+PCOM		TRATE® ROOTT® R - M2A30+PCOM	●	●	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	
	M2A30+PCOM		TRATE® ROOTT® R - M2A45+PCOM	●	●	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	
	M2A45+PCOM		TRATE® ROOTT® R - M3A15+PCOM	●	●	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	
	M3A15+PCOM		TRATE® ROOTT® R - M3A30+PCOM	●	●	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	
	M3A30+PCOM		TRATE® ROOTT® R - M3A45+PCOM	●	●	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	
	M3A45+PCOM		TRATE® ROOTT® R - M4A15+PCOM	●	●	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	
	M4A15+PCOM		TRATE® ROOTT® R - M4A30+PCOM	●	●	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	
	M4A30+PCOM		TRATE® ROOTT® R - M4A45+PCOM	●	●	●	●	●	●	●	●	●	●	●	●	●		●	●	●	●	
	M4A45+PCOM																					